

# **SIMPLYGO MINI Portable Oxygen Concentrator with Extended Battery**

**medsurge**   
*technics for life*



## **SIMPLYGO MINI Portable Oxygen Concentrator with Extended Battery**

### **PRODUCT DESCRIPTION**

One of the most sophisticated and feather-light portable oxygen concentrators on the market today is the SimpleGo MINI by Philips. Weighing in at only 6 pounds with the extended battery attached, it can provide up to 9 hours of therapeutic oxygen.

The SimplyGo MINI delivers pulse flow settings 1 through 5.

## Features

- Attractive and sleek design
- Easy-to-remove external battery
- Light weight – 5.0lbs / With extended Battery is 6 lbs
- Intuitive, easy-to-read screen
- Strong, comfortable carrying case
- Distills air consistently from the environment
- Expunge nitrogen to convey oxygen-rich air
- Provides cleaned oxygen by means of nasal cannula
- Powered by a rechargeable battery or any AC or DC force source.
- Avoid any open flames in the vicinity of oxygen utilization.  
This would incorporate matches, cigarette lighters, candles or smoking in any structure. On the off chance that anybody wishes to smoke while you are utilizing a Portable Oxygen Concentrator, benevolent demand that they venture outside or out of scope of you and the unit.
- Exercise alert around different wellsprings of warmth, for example, compact electric or gas radiators, and additionally broilers and stoves. Keep no less than five feet between such wellsprings of warmth and POC oxygen sources.
- Avoid utilizing petroleum-based moisturizers or creams but instead utilization water-based items.
- Strictly follow the safety and security directions provided on it.
- Always counsel our group of exceedingly expert doctors and health care providers to build up an arrangement that you are alright with or on regions you are not cleared.
- Color Touch Screen LCD
- Powers 4.5 hours on setting 2 and with extended battery 9 hours on setting 2



### **Our SimplyGo MINI with (1) Extended 9 Hour Battery Package Includes:**

- Brand new SimplyGo MINI concentrator
- Oxygen delivery 24/7, daytime and sleep use
- ONE LARGE EXTENDED lithium ion battery (UP TO 9 HOURS)  
(Additional Battery Bundle Sale Packages Available!)
- Convenient carrying bag with comfortable shoulder strap
- Pulse dose oxygen delivery with auto breath detection
- A/C wall power supply- charge your MINI while using
- US and International power supply converters
- D/C vehicle power supply system- Extends Battery Usage
- Easy-read LCD control panel
- Built In "Low O<sub>2</sub>" monitor
- Comfort Plus Soft-tip nasal cannula
- Full manufacturer warranty from Philips  
(Call 1-800-958-0192 for Complete Details)
- Lifetime technical support



### **The SimplyGo MINI comes with two battery options:**

- Smaller standard battery- the MINI unit equipped with the standard battery weighs only 5 pounds and can offer up to 4.5 hours of use at setting 2.
- More powerful extended battery- with this larger battery, the MINI machine weighs in at 6 pounds, and provides up to 9 hours of use at setting 2.

### **Highlights of the SimplyGo MINI with Extended Battery**

- 5 different pulse dose oxygen flow settings
- Lightweight and easy to carry or wear
- Battery life up to 9 hours (double that of the standard battery)
- Easy-to-use touch screen controls
- Hassle-free replacement of rechargeable battery cartridge
- The SimplyGo MINI comes with a stylish carrying case with comfort shoulder strap for your convenience. This oxygen unit is designed to go with you, and maximize your freedom.

### **Available battery options:**

- Small Standard battery – SimplyGo Mini weight of 5 lbs.
- Large Extended battery – SimplyGo Mini weight of 6 lbs.

### **SimplyGo MINI Standard Battery Run Times:**

- Setting 1 = 4.5 Hours
- Setting 2 = 4.0 Hours
- Setting 3 = 2.5 Hours
- Setting 4 = 2.0 Hours
- Setting 5 = 1.5 Hours

### **SimplyGo MINI Extended Battery Run Times:**

- Setting 1 – 9.0 Hours
- Setting 2 – 8.0 Hours
- Setting 3 – 5.0 Hours
- Setting 4 – 4.0 Hours
- Setting 5 – 3.0 Hours

### **Specification:**

- Dimensions: 9.4 x 8.3 x 3.6 Inches with Standard Battery; 10.2 x 8.3 x 2.6 Inches with Extended Battery
- Weight: 5 lbs Standard Battery; 6 lbs Extended Battery
- Control Panel: Color Touch Screen LCD
- Pulse Flow Setting: 1-5
- Universal Power: 100-240 VAC, 50/60 Hz, 2.0 A
- Oxygen Purity/Concentration: 87 to 96% at All Settings
- Low Noise Level: 43 dBA Typical at Setting 2
- Recommended Operating Altitude: Up to 10,000 Feet
- Standard Warranty: 3 Years Device, 1 year Sieve/Accessories



MUSEUM HILL CENTRE, MUSEUM ROAD  
P.O. BOX 75534 00200  
NAIROBI, KENYA  
Tel: 0720 714 337  
Email: [info@medsurgehealth.co.ke](mailto:info@medsurgehealth.co.ke)  
Web: [www.medsurgehealth.co.ke](http://www.medsurgehealth.co.ke)